



Important Infection Prevention

Before beginning this project, please ask yourself the following questions:

Have you been in contact with anyone suspected or confirmed to have COVID-19 in the past 14 days?

Do you have fever, cough or shortness of breath?

If you are able to answer “NO” to both questions, then proceed with making masks.

If you answer “YES” to either or both questions, ***please do not make any masks.***

If you are experiencing cough, fever, or shortness of breath and have a concern that you may have COVID-19 please call your primary care provider for a risk assessment.

Facemask Sewing Instructions

SUPPLIES

100% unused cotton fabric (front) – no metallic fabrics
100% cotton or cotton flannel (back)
1/4” or 3/8” flat elastic

INSTRUCTIONS

Cut the elastic 7” long and tie a knot at each end.
(DO NOT knot the ends of the flat)

Put right sides of cotton fabric together
Cut 9x6 (Adult) or 7.5 x 5 (Child)

Starting at the center of the bottom edge, sew to the first corner, and stop.

Sew the elastic with the edge out into the corner. A few stitches forward and back will hold this.

Sew to the next corner, stop, and bring the other end of the same elastic to the corner and sew a few stitches forward and back.

Now sew across that top of the mask to the next corner. Again put an elastic with the edge out.

Sew to the next corner and sew in the other end of the same elastic.

Sew across the bottom leaving about 1.5” to 2” open. Stop, cut the thread. Turn inside out.

Pin 3 tucks on each side of the mask. Make sure the tucks are the same direction

Sew around the edge of the mask twice.

Masks can be dropped off in re-sealable plastic bags.

For video instructions visit daykimball.org/DonatePPE.
For inquiries or to make arrangements to drop off donations email DonatePPE@daykimball.org.

*Thank you for donating your time and talent.
We appreciate your generosity.*